



Now save  
up to  
**30%\***

## Attractive special prices

Don't miss our special price for  
Transferpette® S micropipette  
from April 1 until June 30, 2024!

**BRAND. For lab. For life.®**



# Promotion – Pipetting

Promotion valid from April 1 until June 30, 2024

## Save up to 30%\*!

Further information from your distributor or at [brand.de/promotions](https://brand.de/promotions).

## Transferpette® S

- + Efficient: one-hand operation
- + Flexible: ergonomic design for all hand sizes
- + Reliable: volume-change protection prevents accidental adjustments

### Transferpette® S, adjustable volume

Items supplied: Transferpette® S, DE-M marking, performance certificate, shelf/rack mount and silicone grease

	Capacity, µl (color-coded)	A** ≤ ± %	µl	CV** ≤ %	µl	Subdivision µl	Cat. No.	Special price GBP*
●	0.1 – 1	2	0.02	1.2	0.012	0.001	705868	185,57
●	0.1 – 2.5	1.4	0.035	0.7	0.018	0.002	705869	185,57
●	0.5 – 10	1	0.1	0.5	0.05	0.01	705870	185,57
●	2 – 20	0.8	0.16	0.4	0.08	0.02	705872	185,57
●	5 – 50	0.8	0.4	0.3	0.15	0.05	705873	185,57
●	10 – 100	0.6	0.6	0.2	0.2	0.1	705874	185,57
●	20 – 200	0.6	1.2	0.2	0.4	0.2	705878	185,57
●	100 – 1000	0.6	6	0.2	2	1	705880	185,57
●	500 – 5000	0.6	30	0.2	10	5	705882	185,57
●	1000 – 10000	0.6	60	0.2	20	10	705884	185,57



### Transferpette® S, fixed volume

Items supplied: Transferpette® S, DE-M marking, performance certificate, shelf / rack mount and silicone grease

	Capacity, µl (color-coded)	A** ≤ ± %	µl	CV** ≤ %	µl	Cat. No.	Special price GBP*
●	10	1	0.1	0.5	0.05	705808	117,85
●	20	0.8	0.16	0.4	0.08	705816	117,85
●	25	0.8	0.2	0.4	0.1	705820	117,85
●	50	0.8	0.4	0.4	0.2	705828	117,85
●	100	0.6	0.6	0.2	0.4	705838	117,85
●	200	0.6	1.2	0.2	0.2	705844	117,85
●	500	0.6	3	0.2	0.2	705854	117,85
●	1000	0.6	6	0.2	0.2	705862	117,85
●	2000	0.8	16	0.3	0.3	705864	117,85



\*\* Calibrated to deliver (TD, Ex). Error limits according to the nominal capacity (= maximum volume) indicated on the instrument, obtained with instrument and distilled water at equilibrium with ambient temperature at 20 °C, and with smooth operation.

A = Accuracy, CV = Coefficient of variation

\* All prices are recommended sales prices without VAT. Savings in percent compared with BRAND price list 2024.

BRAND GMBH + CO KG

P.O. Box 1155 | 97861 Wertheim | Germany

T +49 9342 808 0 | F +49 9342 808 98000 | info@brand.de | www.brand.de



BRAND. For lab. For life.®

BIO-CERT®, BRAND®, BRAND. For lab. For life.®, Transferpette®, as well as the BRAND figurative mark are registered trademarks or trademarks of BRAND GMBH + CO KG, Germany. The BRANDGROUP figurative mark is registered trademark or trademark of Brand Group SE & Co. KG, Germany. All other trademarks mentioned or depicted here are the property of the respective owners.

Subject to technical modification without notice. Errors excepted.

© 2024 BRAND GMBH + CO KG | Printed in Germany | 0424



#### Renewable energy

In our factory we use 100% eco-friendly energy from certified hydro-electric power plants and energy generated on premise from our high efficiency cogeneration unit.



#### Sustainable packaging

We use cardboard with approx 90% recycled content for our product packaging.



#### Waste reduction

The amount of waste is reduced as follows compared to the products mentioned: TipRack compared to TipBox by more than 20%.



BRAND has been awarded **Gold status**

as a sustainable company by **EcoVadis**. For details see [www.brand.de](http://www.brand.de)



A CALIBRE SCIENTIFIC COMPANY