

LABstuff®


Storage

Storage Box 100 & 200	88
MULTI-600 Plastic Tube Racks	89
Stainless Steel Racking	90-91
Stainless Steel Test Tube Racks	92



LabStuff® Storage Box 100

- Box dimensions (includes hinges and front lock): 60 x 138 x 154mm (H x W x D)
- Can accommodate microtubes up to 53.5mm in height.
- 6 colours available - natural, red, green, blue, yellow and black.
- Made from virgin polypropylene, which is autoclavable and translucent, for easy visibility.
- Suitable for 0.5 to 2.0ml volume, capped microtubes and cryovials.
- Grooved top lid for easy stacking.
- Temperature range from -80°C to 135°C.



LabStuff® Storage Box 100

CATALOGUE NUMBER	DESCRIPTION	UNIT QTY	UNIT PRICE (£)
BOX-CRYO-00N	Storage Box 100, natural	5	30.00
BOX-CRYO-00R	Storage Box 100, red	5	30.00
BOX-CRYO-00B	Storage Box 100, blue	5	30.00
BOX-CRYO-00G	Storage Box 100, green	5	30.00
BOX-CRYO-00Y	Storage Box 100, yellow	5	30.00
BOX-CRYO-0BK	Storage Box 100, black	5	30.00

LabStuff® Storage Box 200



- Box dimensions (includes hinges and front lock): 152 x 275 x 62mm
- Polypropylene base and hinged lid.
- Large capacity, 200 place storage box.
- Alphanumeric grid moulded into lid.
- Stackable, in combination with both 100 and 200 place box.
- Autoclavable at 121°C.

LabStuff® Storage Box 200

CATALOGUE NUMBER	DESCRIPTION	UNIT QTY	UNIT PRICE (£)
BOX-CRYO-200	Storage Box 200, natural only	5	60.00

**Stainless Steel
Freezer Racking**

90-91

LabStuff® MULTI-600 Tube Racks



**10%
discount
if you buy
5+ racks**

An attractive range of autoclavable tube racks available in 3 versions.

MULTI-600 is designed to accommodate a full range of laboratory test and centrifuge tubes upto 30mm in diameter.

- Made from highly resistant acetal polymer, which will not stain or shatter in contact with most laboratory chemicals.
- Lightweight, compact and stackable
- A 3-tier design to facilitate the stability and insertion of taller tubes.
- Ideal for Test Tubes, Universals, 15 and 50ml Centrifuge tubes – all in one system.

Up to 13mm Diameter Tubes

CATALOGUE NUMBER	DESCRIPTION	UNIT QTY	UNIT PRICE (£)
175-M600-13B	MULTI-600 Tube Rack, 13mm, 84Pos, Blue	each	14.00
175-M600-13G	MULTI-600 Tube Rack, 13mm, 84Pos, Green	each	14.00
175-M600-13Y	MULTI-600 Tube Rack, 13mm, 84Pos, Yellow	each	14.00

Up to 16mm Diameter Tubes

CATALOGUE NUMBER	DESCRIPTION	UNIT QTY	UNIT PRICE (£)
175-M600-16B	MULTI-600 Tube Rack, 16mm, 60Pos, Blue	each	14.00
175-M600-16G	MULTI-600 Tube Rack, 16mm, 60Pos, Green	each	14.00
175-M600-16Y	MULTI-600 Tube Rack, 16mm, 60Pos, Yellow	each	14.00

Up to 30mm Diameter Tubes

CATALOGUE NUMBER	DESCRIPTION	UNIT QTY	UNIT PRICE (£)
175-M600-30B	MULTI-600 Tube Rack, 30mm, 18Pos, Blue	each	14.00
175-M600-30G	MULTI-600 Tube Rack, 30mm, 18Pos, Green	each	14.00
175-M600-30Y	MULTI-600 Tube Rack, 30mm, 18Pos, Yellow	each	14.00

47

**Racks for
Storage Vials**

LabStuff® SSR - Stainless Steel Racking

Freezer storage solutions

Upright Freezer Racking with Dividers



High-density polypropylene dividers (combs) are used to store virtually any size tube or bottle in the most cost effective and flexible way. Standard division sizes range from 10.5mm to 45mm. Combs are available in 8 different colours to help identify samples.

Trays are offered in either our standard premium grade with welded corners (shown above on left) or cost-saving economy edges (shown on right).

Racks can be 2 to 8 trays high and are available in a variety of sizes to suit most freezers.

Upright Freezer Racking for Boxes



No matter the size box you want to use, we can store it. You can choose from front access with pull-out trays (shown on left) in a variety of formats or a racked, side access system (on right).

Upright Freezer Racking for Microtitre Plates



Open or lidded, shallow or deep, all plates can be accommodated in a made-to-measure rack. Racks can store plates individually or in stacks and can be open fronted, with front retaining rods (to ensure plates remain secure in the rack) or with removable front covers.

LabStuff® SSR - Stainless Steel Racking

Freezer storage solutions

Freezer Bins for General Storage



When you need to store a wide variety of samples, Elkay stainless steel bins can accommodate virtually anything. With several lid options like drop on, hinged or lockable and the added choice of removable or fixed stainless steel dividers, we will help you make the most of your available space.

Silicone-sealed autoclave bins are also available.

Chest Freezer Racking with Dividers

High-density polypropylene dividers (combs) are used to store virtually any size tube or bottle in the most cost effective and flexible way. Combs are available in 8 different colours to help identify samples.

Trays are offered in either our standard premium grade finish with welded corners or cost-saving economy edges (see images at top of page 90). Deep trays are also available for taller tubes.

Front retaining rods are a useful option to ensure trays do not open unintentionally.

Racking is available in full or half height.

Chest Freezer Racking for Boxes and Microtitre Plates



Any size box is catered for as each rack is made-to-measure. Microtitre plates can be stored individually or in stacks.

Racks are available in full or half height, and you also have the choice of: solid or open sides; individually wired compartments or shelves with lipped fronts and front retaining rods to ensure boxes/plates remain secure in the racks.

Petri Dish Racks

Available with open or solid backs, for large or small dishes. Convenient label card holders optimal on the front for easy recognition of samples.



Each of our sales representatives are fully trained in measuring the storage space you have available and can discuss with you all of the options relevant to your needs. Please contact us to arrange a free consultation.

LabStuff® Stainless Steel Test Tube Racks

Freezer storage solutions

Elkay stainless steel test tube racks are ideal when a robust, quality solution for test tube storage is needed, either for use on top of the bench or in a freezer. All racks are autoclavable and are available with two levels, making them ideal for use with round-bottomed tubes and for taller tubes requiring additional support.



Hole Diameter (mm)	14.0	18.0	25.4	30.0	12.0/7.0
Number of Holes	12x5	10x4	6x3	5x2	12x5
Dimensions HxWxD (mm)	268x95x65	268x95x65	268x95x65	268x95x65	268x95x35
Compatible Tube Diameter (mm)	12 & 13	16 & 17	24.5 body 31.0 cap	29.5 body 34.0 cap	1.5ml microcentrifuge
Compatible Tube Height (mm)	75 & 100	100	95	115	n/a
CATALOGUE NUMBER	SSTT-1460	SSTT-1840	SSTT-2518	SSTT-3010	SSTT-1260
Unit Qty	each	each	each	each	each
Unit Price (£)	29.00	29.00	29.00	29.00	29.00

Uses include:
Low-temperature storage, sample preparation, and water baths

Bespoke design service is available to accommodate specific tube or bath sizes