

ExcelStep[®] Repeat Dispenser Tip Selection Guide

Tip Size	Program Number				
	1	2	3	4	5
0,05 ml	1 µl	2 µl	3 µl	4 µl	5 µl
0,50 ml	10 µl	20 µl	30 µl	40 µl	50 µl
1,00 ml	20 µl	40 µl	60 µl	80 µl	100 µl
1,25 ml	25 µl	50 µl	75 µl	100 µl	125 µl
2,50 ml	50 µl	100 µl	150 µl	200 µl	250 µl
5,00 ml	100 µl	200 µl	300 µl	400 µl	500 µl
12,50 ml	250 µl	500 µl	750 µl	1000 µl	1250 µl
25,00 ml	500 µl	1000 µl	1500 µl	2000 µl	2500 µl
50,00 ml	1000 µl	2000 µl	3000 ml	4000 µl	5000 µl
	48	23	15	11	8
	Steps				