

Ripette® Pro Tip Selection Guide

	Program Number									
Tip Size	1	2	3	4	5	6	7	8	9	10
0,05 ml	1 µl	1,5 µl	2 µl	2,5 µl	3 µl	3,5 µl	4 µl	4,5 µl	5 µl	5,5 µl
0,10 ml	2 µl	3 µl	4 µl	5 µl	6 µl	7 µl	8 µl	9 µl	10 µl	11 µl
0,20 ml	4 µl	6 µl	8 µl	10 µl	12 µl	14 µl	16 µl	18 µl	20 µl	22 µl
0,50 ml	10 µl	15 µl	20 µl	25 µl	30 µl	35 µl	40 µl	45 µl	50 µl	55 µl
1,00 ml	20 µl	30 µl	40 µl	50 µl	60 µl	70 µl	80 µl	90 µl	100 µl	110 µl
1,25 ml	25 µl	37,5 µl	50 µl	62,5 µl	75 µl	87,5 µl	100 µl	112,5 µl	125 µl	137,5 µl
2,50 ml	50 µl	75 µl	100 µl	125 µl	150 µl	175 µl	200 µl	225 µl	250 µl	275 µl
5,00 ml	100 µl	150 µl	200 µl	250 µl	300 µl	350 µl	400 µl	450 µl	500 µl	550 µl
10,00 ml	200 µl	300 µl	400 µl	500 µl	600 µl	700 µl	800 µl	900 µl	1000 µl	1100 µl
12,50 ml	250 µl	375 µl	500 µl	625 µl	750 µl	875 µl	1000 µl	1125 µl	1250 µl	1375 µl
25,00 ml	500 µl	750 µl	1000 µl	1250 µl	1500 µl	1750 µl	2000 µl	2250 µl	2500 µl	2750 µl
50,00 ml	1000 µl	1500 µl	2000 ml	2500 µl	3000 µl	3500 µl	4000 µl	4500 µl	5000 µl	5500 µl
	48	35	23	19	15	13	11	9	8	7
	Steps									