


1000 µl (50-1250 µl) OneTouch Barrier - Catalog #10380



Catalog #	Description
10380	Natural, Polypropylene, Sterile

Packaging Type	Packaging Details
Racked Tips	Sterile, 96 Tips/Rack, 8 Racks/Pack, 5 Packs/Case, Hinged Rack

Attributes



Material: Polypropylene
Volume Range: 50-1250 µl
Grad Lines: 200, 500, 1000 µl

Certified RNase/DNase Free
 Certified Human DNA Free
 Certified PCR Inhibitor Free
 Certified Pyrogen Free