


1000 µl (50-1250 µl) OneTouch - Catalog #10360



Catalog #	Description
10360	Natural, Polypropylene
Packaging Type	Packaging Details
Racked Tips	Non-Sterile, 96 Tips/Rack, 8 Racks/Pack, 5 Packs/Case, Hinged Rack
Attributes	



**Volume Range:** 50-1250 µl  
**Grad Lines:** 200, 500, 1000 µl

Certified Human DNA Free  
Certified PCR Inhibitor Free  
Certified Pyrogen Free