


200 µl (1-330 µl) OneTouch Barrier - Catalog #10350



Catalog #

Description

10350

Natural, Polypropylene, Sterile

Packaging Type

Packaging Details

Racked Tips

Sterile, 96 Tips/Rack, 10 Racks/Pack, 5 Packs/Case, Cell Rack

Attributes



Material: Polypropylene

Volume Range: 1-330 µl

Grad Lines: 10, 50, 100 µl

Certified RNase/DNase Free
Certified Human DNA Free
Certified PCR Inhibitor Free
Certified Pyrogen Free