

SciPette

Liquid Handling Solutions

Micro & Multi Pipettors Single/Multichannel



www.sciquip.co.uk

Simply turn the plunger button for volume selection

Finger support keeps the pipettor in place with minimum user effort

Tip ejector allows convenient one-handed operation

Ejector collar and tip cone can be removed for easy cleaning and maintenance

Durable tip cone, piston & ejector collar materials provide excellent chemical resistance



Variable Volume Pipettors

	Part No.	Volume Range	Increment	Test Volume	Inaccuracy (±)	Imprecision (±)
Single-Channel	SP002.5	0.1-2.5µl	0.05µl	2.5µl 1.25µl 0.25µl	2.50% 3.00% 12.00%	2.00% 3.00% 6.00%
	SP0010	0.5-10µl	0.1µl	10µl 5µl 1µl	1.00% 1.50% 2.50%	0.80% 1.50% 1.50%
	SP0020	2-20µl	0.5µl	20µl 10µl 2µl	0.90% 1.20% 3.00%	0.40% 1.00% 2.00%
	SP0050	5-50µl	0.5µl	50µl 25µl 5µl	0.60% 0.90% 2.00%	0.30% 0.60% 2.00%
	SP0100	10-100µl	1µl	100µl 50µl 10µl	0.80% 1.00% 3.00%	0.15% 0.40% 1.50%
	SP0200	20-200µl	1µl	200µl 100µl 20µl	0.60% 0.80% 3.00%	0.15% 0.30% 1.00%
	SP1000A	100-1000µl	5µl	1000µl 800µl 100µl	0.60% 0.70% 2.00%	0.20% 0.25% 0.70%
	SP1000B	200-1000µl	5µl	1000µl 500µl 200µl	0.60% 0.70% 0.90%	0.20% 0.25% 0.30%
	SP5000	1000-5000µl	50µl	5000µl 2500µl 1000µl	0.50% 0.60% 0.70%	0.15% 0.30% 0.30%

SciPette Micro

- Single Channel Pipettors
- Lightweight ergonomic design
- Clear volume display
- Pipettor range from 0.1µl to 5000µl
- Easy to calibrate and maintain
- Designed to reduce physical stress & discomfort
- Serial no. with individual QC certificate
- 3 year warranty

Fixed Volume Pipettors

Part No.	Volume	Test Volume	Inaccuracy (+/-)	Imprecision (+/-)
SP0005FX	5µl	5µl	1.3%	1.2%
SP0010FX	10µl	10µl	0.8%	0.8%
SP0020FX	20µl	20µl	0.6%	0.5%
SP0025FX	250µl	250µl	0.5%	0.3%
SP0050FX	50µl	50µl	0.5%	0.3%
SP0100FX	100µl	100µl	0.5%	0.3%
SP0200FX	200µl	200µl	0.4%	0.2%
SP0250FX	250µl	250µl	0.4%	0.2%
SP0500FX	500µl	500µl	0.3%	0.2%
SP1000FX	1000µl	1000µl	0.3%	0.2%
SP2000FX	2000µl	2000µl	0.3%	0.2%
SP5000FX	5000µl	5000µl	0.3%	0.15%

SciPette Multi

- 8 & 12 channel pipettors
- Ejector works with dispensing head in any position
- Dispensing head rotates for optimum pipetting
- Pipettor range 0.5µl to 300µl
- Individual piston/tip cone assemblies allow easy maintenance
- Compound material tip cone design allows visual seal verification
- Serial no. with individual QC certificate
- 3 year warranty



Multi Channel Pipettors

	Part No.	Volume Range	Increment	Test Volume	Inaccuracy (±)	Imprecision (±)
8-Channel	SP0010-8C	0.5-10µl	0.1µl	10µl 5µl 1µl	1.50% 2.50% 4.00%	1.50% 2.50% 4.00%
	SP0050-8C	5-50µl	0.5µl	50µl 25µl 5µl	1.00% 1.50% 3.00%	0.50% 1.00% 2.00%
	SP0300-8C	50-300µl	5µl	300µl 150µl 50µl	0.70% 1.00% 1.50%	0.25% 0.50% 0.80%
12-Channel	SP0010-12C	0.5-10µl	0.1µl	10µl 5µl 1µl	1.50% 2.50% 4.00%	1.50% 2.50% 4.00%
	SP0050-12V	5-50µl	0.5µl	50µl 25µl 5µl	1.00% 1.50% 3.00%	0.50% 1.00% 2.00%
	SP0300-12C	50-300µl	5µl	300µl 150µl 50µl	0.70% 1.00% 1.50%	0.25% 0.50% 0.80%
Linear Stand	SP2210-LIN	Pipette stand suitable for both single and multi channel pipettors				

SciPette TITAN

- Full volume range 0.1 - 100ml
- Easy to use one-handed operation
- Lightweight, ergonomic design allows longer, fatigue-free pipetting
- Low battery indication
- Fills 100ml quickly with powerful motor



Pipette Filler

Part No.	SP2200
Battery life	15 hrs+ with continuous use
Pipette Types	Glass or Plastic Pipettes (0.1-100ml), Pasteur



SciPette FLOW

- Comfortable and simple to use
- Precise Pipetting control
- Robust and lightweight
- Fits both glass & plastic volumetric & graduated pipettes from 0.1 to 100ml
- Uses replaceable filter

Pipette Controllers

Part No.	Colour
SP2205-BLU	blue
SP2205-RED	red
SP2205-GRE	green
SP2205-YEL	yellow
SP2205-PIN	pink

SciQuip Ltd
Newtown
Wem
Shrewsbury
Shropshire
SY4 5NU, UK

Tel: +44 (0)1939 234222
Fax: +44 (0)1939 234221
Email: info@sciquip.co.uk
www.sciquip.co.uk