

SciPette

Liquid Handling Solutions

Micro & Multi Pipettors
Single/Multichannel



www.sciquip.co.uk

SciPette Multi

- 8 & 12 channel pipettors
- Ejector works with dispensing head in any position
- Dispensing head rotates for optimum pipetting
- Pipettor range 0.5µl to 300µl
- Individual piston/tip cone assemblies allow easy maintenance
- Compound material tip cone design allows visual seal verification
- Serial no. with individual QC certificate
- 3 year warranty



Multi Channel Pipettors

	Part No.	Volume Range	Increment	Test Volume	Inaccuracy (±)	Imprecision (±)
8-Channel	SP0010-8C	0.5-10µl	0.1µl	10µl 5µl 1µl	1.50% 2.50% 4.00%	1.50% 2.50% 4.00%
	SP0050-8C	5-50µl	0.5µl	50µl 25µl 5µl	1.00% 1.50% 3.00%	0.50% 1.00% 2.00%
	SP0300-8C	50-300µl	5µl	300µl 150µl 50µl	0.70% 1.00% 1.50%	0.25% 0.50% 0.80%
12-Channel	SP0010-12C	0.5-10µl	0.1µl	10µl 5µl 1µl	1.50% 2.50% 4.00%	1.50% 2.50% 4.00%
	SP0050-12V	5-50µl	0.5µl	50µl 25µl 5µl	1.00% 1.50% 3.00%	0.50% 1.00% 2.00%
	SP0300-12C	50-300µl	5µl	300µl 150µl 50µl	0.70% 1.00% 1.50%	0.25% 0.50% 0.80%
Linear Stand	SP2210-LIN	Pipette stand suitable for both single and multi channel pipettors				

SciPette TITAN

- Full volume range 0.1 - 100ml
- Easy to use one-handed operation
- Lightweight, ergonomic design allows longer, fatigue-free pipetting
- Low battery indication
- Fills 100ml quickly with powerful motor



Pipette Filler

Part No.	SP2200
Battery life	15 hrs+ with continuous use
Pipette Types	Glass or Plastic Pipettes (0.1-100ml), Pasteur



SciPette FLOW

- Comfortable and simple to use
- Precise Pipetting control
- Robust and lightweight
- Fits both glass & plastic volumetric & graduated pipettes from 0.1 to 100ml
- Uses replaceable filter

Pipette Controllers

Part No.	Colour
SP2205-BLU	blue
SP2205-RED	red
SP2205-GRE	green
SP2205-YEL	yellow
SP2205-PIN	pink

SciQuip Ltd
Newtown
Wem
Shrewsbury
Shropshire
SY4 5NU, UK

Tel: +44 (0)1939 234222
Fax: +44 (0)1939 234221
Email: info@sciquip.co.uk
www.sciquip.co.uk